## 临床领域常用标准菌株目录

中国工业微生物菌种保藏管理中心 孙玉婷 13811635461 sales@china-cicc.org/张京涛 18210192773 zhangjt@china-cicc.org

		icc. org			
Vitek 2 ®Compact	质控菌株	中文名称	MBL货号	Remel货号	CICC编号
	精简质控菌株				
	Enterobacter hormaechei ATCC® 700323 <sup>TM*</sup>	霍氏肠杆菌ATCC® 700323TM*	0755P	R4607094	10432
	Stenotrophomonas maltophilia ATCC® 17666™*	嗜麦芽窄食单胞菌ATCC® 17666TM*	0759P	R4607072	23935
	全面质控菌株				
	Pseudomonas aeruginosa ATCC® 9721 <sup>TM*</sup>	铜绿假单胞菌ATCC® 9721TM*	01009P	R4607068	N/A
	Pseudomonas aeruginosa ATCC® BAA-1744 <sup>TM*</sup>	铜绿假单胞菌ATCC® BAA-1744TM*	01010P	N/A	N/A
GN <b>卡</b>	Acinetobacter baumannii ATCC® BAA-747 <sup>TM*</sup>	鲍曼不动杆菌ATCC® BAA-747™*	0119P	R4601009	N/A
革兰氏阴性细菌鉴定	Ochrobactrum anthropi ATCC® BAA-749™*	人苍白杆菌ATCC® BAA-749™*	0132P	N/A	N/A
	Shigella sonnei ATCC® 25931 <sup>TM*</sup>	宋内志贺氏菌ATCC® 25931TM*	0303P	R4608150	21679
	Proteus vulgaris ATCC® 6380 <sup>TM*</sup>	普通变形菌ATCC® 6380™*	0459P	R4607069	10401
	Enterobacter hormaechei ATCC® 700323™*	霍氏肠杆菌ATCC® 700323TM*	0755P	R4607094	10432
	Klebsiella oxytoca ATCC® 700324 <sup>TM*</sup>	产酸克雷伯菌ATCC® 700324TM*	0757P	R4603912	N/A
	Stenotrophomonas maltophilia ATCC® 17666™*	嗜麦芽窄食单胞菌ATCC® 17666™*	0759P	R4607072	22935
	Elizabethkingia meningoseptica ATCC® 13253 <sup>TM*</sup>	脑膜脓毒性伊丽莎白金菌ATCC® 13253TM*	0971P	R4602002	N/A
	精简质控菌株	•			
	Enterococcus casseliflavus ATCC® 700327 <sup>TM</sup> *	铅黄肠球菌ATCC® 700327TM*	0761P	R4607099	10433
	Staphylococcus saprophyticusATCC® BAA-750 <sup>TM*</sup>	腐生葡萄球菌ATCC® BAA-750™*	0134P	R4607029	N/A
	全面质控菌株	•			
	Streptococcus equi subsp. zooepidemicus ATCC® 43079 <sup>TI</sup> *	□ 马链球菌兽瘟亚种ATCC® 43079™*	0101P	R4608237	N/A
	Kocuria kristinaeATCC® BAA-752 <sup>TM*</sup>	克氏库克菌ATCC® BAA-752™*	0126P	R4603441	N/A
GP卡	Listeria monocytogenesATCC® BAA-751 <sup>TM*</sup>	单核细胞增多李斯特菌ATCC®BAA-751TM*	0130P	R4603971	N/A
革兰氏阳性细菌鉴定	Staphylococcus saprophyticusATCC® BAA-750 <sup>TM*</sup>	腐生葡萄球菌ATCC® BAA-750™*	0134P	R4607029	N/A

-					
	Streptococcus salivarius subsp. thermophilusATCC® 19258 <sup>TM</sup> *	唾液链球菌嗜热亚种ATCC® 19258™ *	0136P	P N/A	6222
	Enterococcus saccharolyticusATCC® 43076 <sup>TM</sup> *	解糖肠球菌ATCC® 43076™*	0223P	N/A R4607099 R4608230 R4609015 R4609015  R4601501 N/A N/A N/A R4601505 R4601501 R4606500 R4601243 R4608012 R4608014 R4604175 R4609374 R4609374 R4601650 N/A R4601602	N/A
	Enterococcus casseliflavusATCC® 700327™*	铅黄肠球菌ATCC® 700327TM*	0761P	R4607099	10433
	Staphylococcus sciuri subsp. sciuriATCC® 29061™*	松鼠葡萄球菌ATCC® 29061TM*	0764P	R4608230	N/A
	Streptococcus pneumoniaeATCC® 49619™*	肺炎链球菌ATCC® 49619™*	0947P	R4609015	N/A
	精简质控菌株	•			
	Candida albicansATCC® 14053™*	白色假丝酵母ATCC® 14053TM*	0332P	R4601501	1944
	全面质控菌株	•			
	Zygosaccharomyces parabailiiATCC® MYA-4549™*	拜氏接合酵母ATCC® MYA-4549™*	01011P	N/A	N/A
	Kloeckera japonica ATCC® 58370™*	日本克勒克酵母ATCC® 58370TM*	01012P	N/A	N/A
	Sporidiobolus salmonicolorATCC® MYA-4550 <sup>TM*</sup>	赭色掷孢酵母ATCC® MYA-4550TM*	01013P	N/A	N/A
YST卡 酵母菌鉴定	Candida glabrataATCC® MYA-2950™*	光滑假丝酵母ATCC® MYA-2950TM*	0122P	R4601505	N/A
	Candida albicansATCC® 14053™*	白色假丝酵母ATCC® 14053TM*	0332P	R4601501	1944
	Staphylococcus epidermidisATCC® 12228 <sup>TM*</sup>	表皮葡萄球菌ATCC® 12228TM*	0371P	R4606500	10294
	Candida lusitaniaeATCC® 34449 <sup>TM*</sup>	葡萄牙假丝酵母ATCC® 34449™*	0774P	R4601243	N/A
	Trichosporon dermatisATCC® 204094™*	粘状丝孢酵母ATCC® 204094TM*	0778P	R4608015	N/A
	Candida utilisATCC® 9950™*	产朊假丝酵母ATCC® 9950TM*	0779P	R4608012	1769
	Prototheca wickerhamiiATCC® 16529™*	魏氏原壁菌ATCC® 16529™*	0780P	R4608014	N/A
	Oligella ureolyticaATCC® 43534™*	解脲寡源杆菌ATCC® 43534TM*	0868P	R4604175	N/A
	精简质控菌株				
	Bacteroides ovatusATCC® BAA-1296™*	卵形拟杆菌ATCC® BAA-1296™*	0585P	R4609374	N/A
	Clostridium septicumATCC® 12464™*	败毒梭菌ATCC® 12464TM*	0586P	R4601650	N/A
	全面质控菌株	•			
	Clostridium perfringensATCC® 13124™*	产气荚膜梭菌ATCC® 13124TM*	0318P	N/A	22949
ANC卡	Clostridium sordelliiATCC® 9714 <sup>TM*</sup>	索氏梭菌ATCC® 9714™*	0331P	R4601602	22950
厌氧菌鉴定	Bacteroides vulgatusATCC® 8482 <sup>TM</sup> *	普通拟杆菌ATCC® 8482™*	0445P	R4601261	10649
	Corynebacterium striatumATCC® BAA-1293™*	纹带棒状杆菌ATCC® BAA-1293™*	0583P	R4609308	N/A

		T	1	1	1
	Parabacteroides distasonisATCC® BAA-1295 <sup>TM*</sup>	迪氏类拟杆菌ATCC® BAA-1295™*	0584P	N/A	N/A
	Bacteroides ovatusATCC® BAA-1296 <sup>TM</sup> *	卵形拟杆菌ATCC® BAA-1296™*	0585P	R4609374	N/A
	Clostridium septicumATCC® 12464 <sup>TM*</sup>	败毒梭菌ATCC® 12464™*	0586P	R4601650	N/A
	精简质控菌株				
	Eikenella corrodensATCC® BAA-1152 <sup>TM*</sup>	啮蚀艾肯菌ATCC® BAA-1152TM*	0189P	R4601947	N/A
	全面质控菌株	-			
	Aggregatibacter aphrophilusATCC® 33389 <sup>TM</sup> *	嗜沫凝聚杆菌ATCC® 33389TM*	0184P	N/A	N/A
*** <del> -</del>	Haemophilus influenzaeATCC® 9007 <sup>TM</sup> *	流感嗜血杆菌ATCC® 9007™*	0185P	R4603815	N/A
VH卡 奈瑟球菌-嗜血杆菌鉴	Eikenella corrodensATCC® BAA-1152 <sup>TM*</sup>	啮蚀艾肯菌ATCC® BAA-1152TM*	0189P	R4601947	N/A
<b></b>	Paenibacillus polymyxaATCC® 7070 <sup>TM</sup> *	多粘类芽孢杆菌ATCC® 7070™*	0296P	R4604600	N/A
	Enterobacter aerogenesATCC® 13048™*	产气肠杆菌ATCC® 13048™*	0306P	R4607080	10293
	Staphylococcus epidermidisATCC® 12228 <sup>TM*</sup>	表皮葡萄球菌ATCC® 12228TM*	0371P	R4606500	10294
	Neisseria gonorrhoeaeATCC® 19424 <sup>TM*</sup>	淋病奈瑟球菌ATCC® 19424TM*	0378P	R4607041	N/A
	Neisseria lactamicaATCC® 23970 <sup>TM</sup> *	乳糖奈瑟球菌ATCC® 23970TM*	0405P	R4607045	N/A
	Oligella urethralisATCC® 17960 <sup>TM*</sup>	尿道寡源杆菌ATCC® 17960TM*	0621P	R4604022	N/A
	Escherichia coliATCC® 25922 <sup>TM*</sup>	大肠埃希菌ATCC® 25922™*	0335P	R4607050	10305
AST-GN卡	Pseudomonas aeruginosaATCC® 27853TM*	铜绿假单胞菌ATCC® 27853TM*	0353P	R4607060	21636
革兰阴性菌药敏试验	Escherichia coliATCC® 35218™*	大肠埃希菌ATCC® 35218TM*	0495P	R4601971	N/A
•	Klebsiella pneumoniae subsp. pneumoniaeATCC® 700603 TM*	肺炎克雷伯菌ATCC® 700603TM*	0784P	R4603074	N/A
	Enterococcus faecalisATCC® 29212 <sup>TM*</sup>	粪肠球菌ATCC® 29212™*	0366P	R4607030	23658
	Enterococcus faecalisATCC® 51299 <sup>TM*</sup>	粪肠球菌ATCC® 51299™*	0959P	R4601996	24243
A CET. CED. H	Staphylococcus aureus subsp. aureusATCC® 29213 <sup>TM*</sup>	金黄色葡萄球菌ATCC® 29213TM*	0365P	R4607011	10473
AST-GP卡 革兰阳性菌药敏试验	Staphylococcus aureus subsp. aureusATCC® BAA-976™*	金黄色葡萄球菌ATCC® BAA-976™*	0146P	R4606512	N/A
	Staphylococcus aureus subsp. aureusATCC® BAA-977 <sup>TM*</sup>	金黄色葡萄球菌ATCC® BAA-977™*	0147P	R4606513	N/A
	Staphylococcus aureusATCC® BAA-1026 <sup>TM</sup> *	金黄色葡萄球菌ATCC® BAA-1026™*	0179P	R4601026	N/A
	Streptococcus pneumoniaeATCC® 49619 <sup>TM</sup> *	肺炎链球菌ATCC® 49619™*	0947P	R4609015	N/A
AST-ST卡 连球菌药敏试验	Staphylococcus aureus subsp. aureusATCC® 29213 <sup>TM*</sup>	金黄色葡萄球菌ATCC® 29213TM*	0365P	R4607011	10473

	Staphylococcus aureus subsp. aureusATCC® BAA-977 <sup>TM</sup> *	金黄色葡萄球菌ATCC® BAA-977™*	0147P	R4606513	N/A
AST-YS卡	Issatchenkia orientalisATCC® 6258 <sup>TM</sup> *	克柔假丝酵母ATCC® 6258™*	0227P	R4601520	N/A
酵母菌药敏试验	Candida parapsilosisATCC® 22019™*	近平滑假丝酵母ATCC® 22019TM*	0726P	R4601518	1973
Vitek®MS	质控菌株	中文名称	MBL货号	Remel货号	CICC编号
细菌	Enterobacter aerogenes ATCC® 13048 <sup>TM</sup> *	产气肠杆菌ATCC® 13048™*	0306P	R4607080	10293
酵母菌	Candida glabrata ATCC® MYA-2950 <sup>TM*</sup>	光滑假丝酵母ATCC® MYA-2950™*	0122P	R4601505	N/A
Mycobacterium和 Nocardia	Mycobacterium smegmatis ATCC® 19420™*	耻垢分枝杆菌ATCC® 19420™*	0514P	N/A	N/A
	Nocardia farcinica Trevisan (ATCC® 3308 <sup>TM</sup> )	皮疽诺卡氏菌ATCC® 3308™	N/A	R4607156	N/A
API生化试验条	质控菌株	中文名称	MBL货号	Remel货号	CICC编号
API® 20E	Enterobacter cloacae subsp. cloacae ATCC® 13047™*	阴沟肠杆菌阴沟亚种ATCC® 13047™	0323P	R4607090	21539
	Escherichia coli ATCC® 25922 <sup>TM</sup> *	大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
	Stenotrophomonas maltophilia ATCC® 51331 <sup>TM</sup> *	嗜麦芽寡养单胞菌ATCC® 51331™*	0742P	R4607064	N/A
	Klebsiella pneumoniae subsp. pneumoniae ATCC® 35657 TM*	肺炎克雷伯氏菌肺炎亚种ATCC® 35657 <sup>TM</sup> *	0942P	R4603922	N/A
	Proteus mirabilis ATCC® 35659 <sup>TM</sup> *	奇异变形菌ATCC® 35659™*	0944P	R4607077	N/A
	Proteus hauseri ATCC® 13315 <sup>TM*</sup>	豪氏变形菌ATCC® 13315™*	0355P	R4607057	10310
API® Rapid 20E	Escherichia coli ATCC® 11775™*	大肠埃希菌ATCC® 11775TM*	0465P	R4601989	24085
•	Klebsiella pneumoniae subsp. pneumoniae ATCC® 35657 TM*	肺炎克雷伯氏菌肺炎亚种ATCC® 35657™*	0942P	R4603922	N/A
	Pseudomonas aeruginosa ATCC® 27853 <sup>TM*</sup>	铜绿假单胞菌ATCC® 27853TM*	0353P	R4607060	21636
	Aeromonas hydrophila ATCC® 35654 <sup>TM</sup> *	嗜水气单胞菌ATCC® 35654TM*	0910P	R4601013	N/A
API® 20 NE	Alcaligenes faecalis subsp. faecalis ATCC® 35655TM*	粪产碱菌ATCC® 35655TM*	0911P	R4601011	10289
	Sphingobacterium multivorum ATCC® 35656™*	多食鞘氨醇杆菌ATCC® 35656TM*	0948P	R4602004	N/A
	Aggregatibacter aphrophilus ATCC® 49917 <sup>TM*</sup>	嗜泡团聚杆菌ATCC® 49917TM*	0249P	R4603902	N/A
API® NH	Neisseria gonorrhoeae ATCC® 31426 <sup>TM</sup> *	淋病奈瑟氏菌ATCC® 31426TM*	0375P	R4607042	N/A
	Haemophilus influenzae (type b) ATCC® 10211™*	流感嗜血杆菌ATCC® 10211TM*	0441P	R4603810	N/A
	Staphylococcus lentus ATCC® 700403 <sup>TM*</sup>	慢葡萄球菌ATCC® 700403TM*	0739P	R4606511	N/A
API®Staph	Staphylococcus xylosus ATCC® 700404TM*	木糖葡萄球菌ATCC® 700404™*	0741P	R4606506	N/A

	Staphylococcus capitis subsp. capitis ATCC® 35661 <sup>TM*</sup>	头葡萄球菌头亚种ATCC® 35661TM*	0963P	R4608207	10290
	T ,	>			
RAPIDEC® Staph	Staphylococcus aureus subsp. aureus ATCC® 25923™*	金黄色葡萄球菌ATCC® 25923TM*	0360P	R4607010	23926
ин и изе с вири	Staphylococcus saprophyticus ATCC® 43867™*	腐生葡萄球菌ATCC® 43867™*	N/A	R4608223	N/A
API® 20 Strep	Streptococcus equi subsp. zooepidemicus ATCC® 700400 тм*	马链球菌兽瘟亚种ATCC® 700400™*	0743P	R4608225	N/A
	Streptococcus uberis ATCC® 700407™*	乳房链球菌ATCC® 700407™*	0744P	R4608236	N/A
	Corynebacterium renale ATCC® 19412 <sup>TM*</sup>	牛肾盂炎棒杆菌ATCC® 19412TM*	0293P	R4601906	N/A
.PI®Coryne	Cellulosimicrobium cellulans ATCC® 27402™*	纤维纤维微菌ATCC® 27402™*	0294P	R4601544	N/A
	Microbacterium testaceum ATCC® 15829™*	砖红色微杆菌ATCC® 15829™*	0295P	R4604040	N/A
	Clostridium perfringens ATCC® 13124™*	产气荚膜梭菌ATCC® 13124™*	0318P	N/A	22949
API® 20 A	Clostridium sordellii ATCC® 9714™*	索氏梭菌ATCC® 9714™*	0331P	R4601602	22950
	Bacteroides ovatus ATCC® 8483™*	卵形拟杆菌ATCC® 8483™*	0400P	R4601222	N/A
	Clostridium sporogenes ATCC® 19404™*	生孢梭菌ATCC® 19404™*	0317P	R4601700	10409
	Clostridium sordellii ATCC® 9714™*	索氏梭菌ATCC® 9714™*	0331P	R4601602	22950
Rapid ID32A	Capnocytophaga sputigena ATCC® 33612™*	生痰二氧化碳嗜纤维菌ATCC® 33612 TM*	0749P	N/A	N/A
	Actinomyces viscosus ATCC® 15987 <sup>TM*</sup>	黏放线菌ATCC® 15987TM*	0750P	N/A	N/A
	Bacteroides fragilis ATCC® 23745 <sup>TM*</sup>	脆弱拟杆菌ATCC® 23745TM*	0940P	R4601251	N/A
	Cryptococcus laurentii ATCC® 18803™*	劳伦梯氏隐球菌ATCC® 18803™*	0333P	R4601241	N/A
API® 20 C AUX	Candida glabrata ATCC® 15126 <sup>TM*</sup>	光滑假丝酵母ATCC® 15126™*	0737P	R4601513	N/A
	Meyerozyma guilliermondii ATCC® 6260™*	季也蒙毕赤酵母ATCC® 6260™*	0738P	R4601521	1945
API®ZYM	Pseudomonas aeruginosa ATCC 27853	铜绿假单胞菌ATCC 27853	0353P	R4607060	21636
	Paenibacillus polymyxa ATCC® 43865™*	多粘类芽孢杆菌ATCC® 43865™*	0228P	R4601234	10288
API®50CH	Lactobacillus paracasei subsp. paracasei ATCC® BAA-52 TM*	类干酪乳杆菌类干酪亚种ATCC® BAA-52™*	0526P	N/A	N/A
	Klebsiella pneumoniae subsp. pneumoniae ATCC® 35657 TM*	肺炎克雷伯氏菌肺炎亚种ATCC® 35657™*	0942P	R4603922	N/A
霉菌	Aspergillus brasiliensis ATCC® 16404 <sup>TM*</sup>	巴西曲霉ATCC® 16404™*	0392P	R4601100	2463
培养基种类	质控菌株	中文名称	MBL货号	Remel货号	CICC编号

营养琼脂培养基	Staphylococcus aureus subsp. aureus ATCC® 25923 <sup>TM</sup> *	金黄色葡萄球菌ATCC® 25923TM*	0360P	R4607010	23926
	Pseudomonas aeruginosa ATCC® 27853 <sup>TM</sup> *	## ATCC® 27853TM*	0353P	R4607060	21636
血琼脂培养基	Streptococcus pyogenes ATCC® 19615 <sup>TM*</sup>	化脓链球菌ATCC® 19615™*	0385P	R4607000	10464
皿坏加切介坐	Streptococcus pneumoniae ATCC® 6305TM*	肺炎链球菌ATCC® 6305™*	0447P	R4607024	N/A
	Neisseria gonorrhoeae ATCC® 43069 <sup>TM*</sup>	淋病奈瑟氏菌ATCC® 43069TM*	0429P	R4607043	N/A
巧克力培养基	Haemophilus influenzae (type b) ATCC® 10211 <sup>TM</sup> *	流感嗜血杆菌ATCC® 10211TM*	0441P	R4603810	N/A
淋球菌培养基	Neisseria gonorrhoeae derived from ATCC® 49226 <sup>TM*</sup>	淋病奈瑟氏菌ATCC® 49226TM*	0648P	R4609006	N/A
	Escherichia coli ATCC® 25922 <sup>TM*</sup>	大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
M-H培养基	Staphylococcus aureus subsp. aureus ATCC® 25923 <sup>TM</sup> *	金黄色葡萄球菌ATCC® 25923TM*	0360P	R4607010	23926
	Escherichia coli ATCC® 25922 <sup>TM*</sup>	大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
HTM培养基	Staphylococcus aureus subsp. aureus ATCC® 25923™*	金黄色葡萄球菌ATCC® 25923TM*	0360P	R4607010	23926
营养肉汤培养基	Escherichia coli ATCC® 25922™*	大肠埃希菌ATCC® 25922™*	0335P	R4607050	10305
	Staphylococcus aureus subsp. aureus ATCC® 25923™*	金黄色葡萄球菌ATCC® 25923TM*	0360P	R4607010	23926
牛脑心浸液培养基	Escherichia coli ATCC® 25922™*	大肠埃希菌ATCC® 25922™*	0335P	R4607050	10305
	Staphylococcus aureus subsp. aureus ATCC® 25923™*	金黄色葡萄球菌ATCC® 25923TM*	0360P	R4607010	23926
	Bacteroides fragilis ATCC® 25285TM*	脆弱拟杆菌ATCC® 25285™*	0320P	R4601250	24309
厌氧血琼脂培养基	Prevotella melaninogenica ATCC® 25845 <sup>TM</sup> *	产黑普氏菌ATCC® 25845TM*	0110P	R4601255	N/A
罗氏培养基	Mycobacterium tuberculosisATCC® 25177 <sup>TM*</sup>	结核分枝杆菌ATCC® 25177TM*	0112P	N/A	N/A
	Mycobacterium kansasii ATCC® 12478 <sup>TM*</sup>	堪萨斯分枝杆菌ATCC® 12478TM*	0545P	N/A	N/A
染色方法	质控菌株	中文名称	MBL货号	Remel货号	CICC编号
革兰染色	阳性: Staphylococcus aureus subsp. aureus ATCC® 25923™*		0360P	R4607010	23926
半二米口	阴性: Escherichia coli ATCC® 25922TM*	阴性: 大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
抗酸染色	阳性: Mycobacterium chelonae CMCC 93326	阳性: 龟分枝杆菌 CMCC 93326	N/A	N/A	N/A
加政策也	阴性: Escherichia coli ATCC® 25922™*	阴性: 大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
<b>鞴毛</b> 沈	阳性: Pseudomonas aeruginosa ATCC® 27853™*	阳性:铜绿假单胞菌ATCC® 27853TM *	0353P	R4607060	21636

牧七米巳	阴性: Acinetobacter baumannii ATCC® 19606™*	阴性: 鲍氏不动杆菌ATCC® 19606™	0357P	R4601007	N/A
	1911: Acinctobacter baumanin ATCC 19000	*	03371	K4001007	IV/A
墨汁染色	阳性: Cryptococcus neoformans ATCC® 14116 <sup>TM*</sup>	阳性:新生隐球菌 ATCC® 14116™	0291P	R4601812	N/A
	阴性:Candida albicans ATCC® 90028™*	阴性: 白色念珠菌 ATCC® 90028TM*	0264P	R4601496	N/A
试验方法	质控菌株	中文名称	MBL货号	Remel货号	CICC编号
触酶试验	阳性: Staphylococcus aureus subsp. aureus ATCC® 25923™*	阳性: 金黄色葡萄球菌ATCC® 25923 TM*	0360P	R4607010	23926
	阴性: Streptococcus pneumoniaeATCC® 49619™*	阴性: 肺炎链球菌ATCC® 49619™*	0947P	R4609015	N/A
复 / L 耐油	阳性: Pseudomonas aeruginosa ATCC® 27853TM*	阳性:铜绿假单胞菌ATCC® 27853TM	0353P	R4607060	21636
氧化酶试验	阴性: Escherichia coli ATCC® 25922™*	阴性: 大肠埃希菌ATCC® 25922™*	0335P	R4607050	10305
	阳性: Staphylococcus aureus subsp. aureus ATCC®	阳性:金黄色葡萄球菌ATCC® 25923	0360P	R4607010	23926
凝固酶试验	阴性: Staphylococcus epidermidis ATCC® 51625™*	阴性: 表皮葡萄球菌ATCC® 51625™ *	01068P	N/A	N/A
DNA酶试验	阳性: Serratia marcescens ATCC® 14041™*	阳性: 黏质沙雷氏菌ATCC® 14041™	0506P	N/A	10355
DNA時 风驰	阴性: Escherichia coli ATCC® 25922™*	阴性: 大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
CAMP试验	阳性: Streptococcus agalactiae ATCC® 13813™*	阳性: 无乳链球菌ATCC® 13813TM*	0370P	R4608250	10465
CAMPTA 3 <u>w</u>	阴性: Streptococcus pyogenes ATCC® 19615™*	阴性: 化脓链球菌ATCC® 19615™*	0385P	R4607000	10464
β -内酰胺酶试验	阳性: Staphylococcus aureus subsp. aureus ATCC®	阳性:金黄色葡萄球菌ATCC® 29213	0365P	R4607011	10473
<b>户</b> - 内肌放射 风湿	阴性: Staphylococcus aureus subsp. aureus ATCC®	阴性:金黄色葡萄球菌ATCC® 25923	0360P	R4607010	23926
Optochin试验	阳性: Streptococcus pneumoniae ATCC® 49619™*	阳性: 肺炎链球菌ATCC® 49619™*	0947P	R4609015	N/A
Optociiii क्रिज्	阴性: Enterococcus faecalis ATCC® 29212™	阴性: 粪肠球菌ATCC® 29212™*	0366P	R4607030	23658
嗜血杆菌因子试验	阳性: Haemophilus influenzae ATCC® 49247™*	阳性: 流感嗜血杆菌ATCC® 49247™	0647P	R4603830	N/A
增皿杆图囚丁砜巡	阴性: Aggregatibacter aphrophilus ATCC® 7901™*	阴性:嗜泡团聚杆菌ATCC® 7901™*	0411P	R4604000	N/A
糖 (醇、苷) 类发酵	阳性: Escherichia coli ATCC® 25922™*	阳性: 大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
试验	阴性: Morganella morganii subsp. morganii ATCC®	阴性:摩氏摩根氏菌摩根亚种ATCC®	0215P	N/A	N/A
葡萄糖O/F试验	Pseudomonas aeruginosa ATCC® 27853 <sup>TM*</sup>	铜绿假单胞菌ATCC® 27853™*	0353P	R4607060	21636
用街/店U/F以沙	Escherichia coli ATCC® 25922™*	大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
	Escherichia coli ATCC® 25922™*	大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
克氏双糖铁试验	Salmonella enterica subsp. enterica serovar Typhimurium ATCC® 14028 <sup>TM*</sup>	肠道沙门氏菌肠道亚种鼠伤寒血清型ATCC® 14028™*	0363P	R4606000	22956
三糖铁试验	Escherichia coli ATCC® 25922 <sup>TM*</sup>	大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
二個钦阺娅	Shigella flexneri CMCC 51573	弗氏志贺氏菌CMCC 51573	N/A	N/A	N/A
田甘仁汁心	阳性: Escherichia coli ATCC® 25922™*	阳性: 大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
甲基红试验	阴性: Enterobacter cloacae subsp. cloacae ATCC® 13047	阴性: 阴沟肠杆菌阴沟亚种ATCC®	0323P	R4607090	21539

VP试验	阳性: Klebsiella pneumoniae subsp. pneumoniae ATCC® 13883 <sup>TM*</sup>	阳性: 肺炎克雷伯氏菌肺炎亚种ATCC ® 13883 <sup>TM*</sup>	0351P	R4607037	24165
	阴性: Escherichia coli ATCC® 25922™*	阴性: 大肠埃希菌ATCC® 25922™*	0335P	R4607050	10305
1-11 <del>**</del>	阳性: Enterococcus faecalis ATCC® 29212™*	阳性: 粪肠球菌ATCC® 29212TM*	0366P	R4607030	23658
七叶苷试验	阴性: Streptococcus pneumoniae ATCC® 49619™*	阴性: 肺炎链球菌ATCC® 49619™*	0947P	R4609015	N/A
nalinn. v-PriA	阳性: Escherichia coli ATCC® 25922TM*	阳性: 大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
吲哚试验	阴性: Klebsiella pneumoniae subsp. pneumoniae ATCC®	阴性: 肺炎克雷伯氏菌肺炎亚种ATCC	0351P	R4607037	24165
消酸盐还原实验	阳性: Staphylococcus aureus subsp. aureus ATCC®	阳性:金黄色葡萄球菌ATCC® 25923	0360P	R4607010	23926
	阴性: Acinetobacter baumannii ATCC® 19606™*	阴性: 鲍氏不动杆菌ATCC® 19606™	0357P	R4601007	N/A
	阳性: Proteus vulgaris ATCC® 8427™*	阳性: 普通变形菌ATCC® 8427™*	0691P	R4607058	N/A
1水11年1八5至	阴性: Escherichia coli ATCC® 25922™*	阴性: 大肠埃希菌ATCC® 25922™*	0335P	R4607050	10305
干菌肽敏感试验	阳性: Streptococcus pyogenes ATCC® 19615™*	阳性: 化脓链球菌ATCC® 19615™*	0385P	R4607000	10464
竹囷枞蚁怒试验	阴性: Streptococcus agalactiae ATCC® 12386™*	阴性: 无乳链球菌ATCC® 12386™*	0439P	R4607027	N/A
	Shigella sonnei ATCC® 25931 <sup>TM*</sup>	宋内志贺氏菌ATCC® 25931TM*	0303P	R4608150	21679
土加弗布建兴城市	Shigella boydii ATCC® 9207 <sup>TM</sup> *	鲍氏志贺氏菌ATCC® 9207™*	0349P	R4608110	21680
志贺菌血清学鉴定	Shigella flexneri ATCC® 12022 <sup>TM*</sup>	弗氏志贺氏菌ATCC® 12022TM*	0356P	R4608101	21678
	Shigella dysenteriae ATCC® 13313TM*	痢疾志贺氏菌ATCC® 13313TM*	01235P	R4659100	N/A
	Salmonella enterica subsp. enterica serovar Paratyphi A	甲型副伤寒沙门菌A ATCC® 9150™*	0341P	N/A	N/A
	Salmonella enterica subsp. enterica serovar Tallahassee ATCC® 12002 <sup>TM</sup> *	肠道沙门氏菌肠道亚种塔拉哈西血清型ATCC® 12002™*	0344P	N/A	N/A
沙门菌血清学鉴定	Salmonella enterica subsp. enterica serovar Anatum ATCC ® 9270 <sup>TM</sup> *	肠道沙门氏菌肠道亚种鸭血清型ATCC® 9270 <sup>TM*</sup>	0346P	R4606006	N/A
	Salmonella enterica subsp. enterica serovar Typhimurium ATCC® 14028 <sup>TM*</sup>	肠道沙门氏菌肠道亚种鼠伤寒血清型ATCC® 14028™*	0363P	R4606000	22956
	Salmonella enterica subsp. enterica serovar Choleraesuis ATCC® 10708 <sup>TM</sup> *	肠道沙门氏菌肠道亚种猪霍乱血清型 ATCC® 10708™*	0902P	R4609001	N/A
致病性大肠杆菌血清 学鉴定	阳性: 致病性大肠埃希氏菌 Escherichia coli EHEC O157:H7	出血性大肠埃希氏菌CICC 21530	N/A	N/A	21530
<b>子金</b> 化	阴性: Escherichia coli ATCC® 25922™*	大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
CLSI M100	质控菌株	中文名称	MBL货号	Remel货号	CICC编号
	Escherichia coli ATCC® 25922 <sup>TM*</sup>	大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
	Staphylococcus aureus ATCC® 25923TM*	金黄色葡萄球菌ATCC® 25923TM*	0360P	R4607010	23926
非苛养微生物纸片扩 散法	Pseudomonas aeruginosa ATCC® 27853™*	铜绿假单胞菌ATCC® 27853TM*	0353P	R4607060	21636
IXIA	Escherichia coli ATCC® 35218 <sup>TM</sup> *	大肠埃希菌ATCC® 35218TM*	0495P	R4601971	N/A
	Klebsiella pneumoniae ATCC® 700603™*	肺炎克雷伯菌ATCC® 700603™*	0784P	R4603074	N/A
	Haemophilus influenzae ATCC® 49247 <sup>TM</sup> *	流感嗜血杆菌ATCC® 49247™*	0647P	R4603830	N/A

	Neisseria gonorrhoeae ATCC® 49226 <sup>TM</sup> *	淋病奈瑟氏菌ATCC® 49226TM*	0648P	R4609006	N/A
苛养微生物纸片扩散 法	Haemophilus influenzae ATCC® 49766 <sup>TM</sup> *	流感嗜血杆菌ATCC® 49766TM*	0919P	R4603806	N/A
非苛养微生物的MIC测定 持养微生物的MIC测定 林病奈瑟菌的MIC测 医氧菌MIC的测定 (琼脂稀释法和微量 对汤稀释法) 超广谱β内酰胺酶 ESBL)的检测 效良Hodge实验/Carba	Streptococcus pneumoniae ATCC® 49619 <sup>TM</sup> *	肺炎链球菌ATCC® 49619™*	0947P	R4609015	N/A
	Staphylococcus aureus ATCC® 29213 <sup>TM*</sup>	金黄色葡萄球菌ATCC® 29213TM*	0365P	R4607011	10473
	Enterococcus faecalis ATCC® 29212TM*	粪肠球菌ATCC® 29212TM*	0366P	R4607030	23658
非苛养微生物的MIC	Escherichia coli ATCC® 25922™*	大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
测定	Pseudomonas aeruginosa ATCC® 27853™*	铜绿假单胞菌ATCC® 27853TM*	0353P	R4607060	21636
	Escherichia coli ATCC® 35218 <sup>TM</sup> *	大肠埃希菌ATCC® 35218TM*	0495P	R4601971	N/A
	Klebsiella pneumoniae ATCC® 700603™*	肺炎克雷伯菌ATCC® 700603™*	0784P	R4603074	N/A
<b></b> 5 5 5 8 8 9 8 9 9 9 9 9 9 9 9 9 9 9 9 9	Haemophilus influenzae ATCC® 49247™*	流感嗜血杆菌ATCC® 49247™*	0647P	R4603830	N/A
苛养微生物的MIC测定 定 淋病奈瑟菌的MIC测	Haemophilus influenzae ATCC® 49766™*	流感嗜血杆菌ATCC® 49766™*	0919P	R4603806	N/A
	Streptococcus pneumoniae ATCC® 49619 <sup>TM</sup> *	肺炎链球菌ATCC® 49619™*	0947P	R4609015	N/A
淋病奈瑟菌的MIC测 定	Neisseria gonorrhoeae derived from ATCC® 49226™*	淋病奈瑟氏菌ATCC® 49226™*	0648P	R4609006	N/A
	Bacteroides fragilis ATCC® 25285 <sup>TM*</sup>	脆弱拟杆菌ATCC® 25285TM*	0320P	R4601250	24309
厌氧菌MIC的测定 (克斯稀释注和微量	Bacteroides thetaiotaomicron ATCC® 29741™*	多形拟杆菌ATCC® 29741™*	0319P	R4601260	N/A
肉汤稀释法)	Clostridium difficile ATCC® 700057TM*	艰难梭菌ATCC® 700057™*	0527P	N/A	N/A
	Eggerthella lenta ATCC® 43055TM*	缓慢艾格士氏菌ATCC® 43055TM*	0936P	R4601951	N/A
起亡途Q 力融 腔 施	Escherichia coli ATCC® 25922™*	大肠埃希菌ATCC® 25922TM*	0335P	R4607050	10305
(ESBL)的检测	Klebsiella pneumoniae subsp. pneumoniae ATCC® 700603 TM*	肺炎克雷伯菌ATCC® 700603™*	0784P	R4603074	N/A
改良Hodge实验/Carba	Klebsiella pneumoniae ATCC® BAA-1705 <sup>TM*</sup>	肺炎克雷伯菌ATCC® BAA-1705™*	01005P	N/A	N/A
NP试验检测碳青霉烯 酶	Klebsiella pneumoniae ATCC® BAA-1706™*	肺炎克雷伯菌ATCC® BAA-1706™*	01006P	N/A	N/A
头孢硝噻吩纸片法检	Staphylococcus aureus subsp. aureus ATCC® 25923™*	金黄色葡萄球菌ATCC® 25923TM*	0360P	R4607010	23926
则β -内酰胺酶	Staphylococcus aureus subsp. aureus ATCC® 29213 <sup>TM*</sup>	金黄色葡萄球菌ATCC® 29213TM*	0365P	R4607011	10473
耐甲氧西林金黄色葡 萄球菌(MRSA)检	Staphylococcus aureus subsp. aureus ATCC® 29213 <sup>TM*</sup>	金黄色葡萄球菌ATCC® 29213TM*	0365P	R4607011	10473
测	Staphylococcus aureus subsp. aureus ATCC® 43300™*	金黄色葡萄球菌ATCC® 43300™*	0852P	R4609022	N/A

BD Phoenix	质控菌株	中文名称	MBL货号	Remel货号	CICC编号
NMIC/ID药敏+鉴定	Escherichia coliATCC® 25922 <sup>TM*</sup>	大肠埃希菌ATCC® 25922™*	0335P	R4607050	10305
	Pseudomonas aeruginosa ATCC® 27853™*	铜绿假单胞菌ATCC® 27853TM*	0353P	R4607060	21636
	Escherichia coli ATCC® 35218 <sup>TM</sup> *	大肠埃希菌ATCC® 35218TM*	0495P	R4601971	N/A
	Klebsiella pneumoniae subsp. pneumoniae ATCC® 700603 TM*	肺炎克雷伯菌ATCC® 700603TM*	0784P	R4603074	N/A
	Staphylococcus aureus subsp. aureus ATCC® 25923 <sup>TM*</sup>	金黄色葡萄球菌ATCC® 25923TM*	0360P	R4607010	23926
	Staphylococcus aureus subsp. AureusATCC® 29213 <sup>TM*</sup>	金黄色葡萄球菌ATCC® 29213TM*	0365P	R4607011	10473
PMIC/ID药敏+鉴定	Staphylococcus aureus ATCC® BAA-976™*	金黄色葡萄球菌ATCC® BAA-976TM*	0146P	R4606512	N/A
	Staphylococcus aureus ATCC® BAA-977 <sup>TM</sup> *	金黄色葡萄球菌ATCC® BAA-977™*	0147P	R4606513	N/A
	Enterococcus faecalis ATCC® 29212 <sup>TM</sup> *	粪肠球菌ATCC® 29212TM*	0366P	R4607030	23658
	Enterococcus faecalis ATCC® 51299 <sup>TM</sup> *	粪肠球菌ATCC® 51299™*	0959P	R4601996	24243
SMIC/ID药敏+鉴定	Streptococcus agalactiae ATCC® 13813TM*	无乳链球菌ATCC® 13813™*	0370P	R4608250	10465
	Streptococcus pneumoniae ATCC® 49619 <sup>TM*</sup>	肺炎链球菌ATCC® 49619™*	0947P	R4609015	N/A